



ADFLEX- I.V 에드플렉스- I.V

(PU WATERPOOF DRESSING)



INTENDED USE

Used to fix the needle, fix the device placed under the skin or in the blood vessel.




SPECIFICATIONS

- Polyurethane (P.U) dressing
- Cover with Acrylic glue
- Stickiness 440g/12mm $\pm$ 30g, Tensile strength 2.0kg/12mm $\pm$ 200g
- Quality standards/ Certificate: ISO 13485, CE, FDA.

PRODUCT CHARACTERISTICS

- **Polyurethane (P.U) tape**, Latex-free, **transparent**, easy to monitor the tape area, waterproof.
- **Breathable**, semi-permeable P.U film, allowing oxygen exchange and moisture loss, but not allowing bacteria and water to pass through (waterproof) to help protect the tape area from environmental impacts, 4 rounded corners to prevent edges from coming loose.
- **There is a special groove to fix the needle.**
- **Acrylic glue** has good adhesion and **is safe for all skin types.**
- **Each piece is sterilized with Ethylene Oxide (E.O Gas).**
- **Sterile packaging for each piece individually.**

HOW TO USE



ITEM	SIZE	GROOVE	PACKING
AF5160	5.1 cm x 6cm	4mm x 25mm	100 pcs/ Box
AFI6710	6.7 cm x 10cm (oval)	7mm x 25mm	50 pcs/ Box

Băng keo cố định kim luôn
không thấm nước



ADFLEX- I.V

(P.U WATERPOOF DRESSING)



CHỈ ĐỊNH

Sử dụng để cố định kim luôn, cố định thiết bị đặt dưới da hoặc trong mạch máu.




THÔNG SỐ

- Băng Polyurethane (P.U)
- Phủ keo Acrylic
- Độ dính 440g/12mm $\pm$ 30g, Độ bền kéo 2.0kg/12 mm $\pm$ 200g
- Tiêu chuẩn chất lượng/ Giấy chứng nhận: ISO 13485, CE, FDA.

ĐẶC TÍNH SẢN PHẨM

- Băng Polyurethane (P.U), không chứa Latex, trong suốt tạo điều kiện dễ dàng theo dõi vùng dán băng, không thấm nước.
- Màng film P.U thoáng khí, bán thấm, cho phép trao đổi Oxy thoát hơi ẩm, nhưng không cho vi khuẩn và nước đi qua (không thấm nước) giúp bảo vệ vùng dán băng khỏi tác động của môi trường, 4 góc bo tròn tránh bung mép.
- Có đường xé rãnh chuyên dùng cố định kim luôn.
- Keo Acrylic độ dính tốt và an toàn cho mọi loại da.
- Tiệt trùng từng miếng bằng ethylene oxide (E.O Gas).
- Đóng gói vô trùng riêng từng miếng.

HƯỚNG DẪN SỬ DỤNG



Mã sản phẩm	Kích cỡ	Đường xé rãnh	Đóng gói
AF5160	5.1 cm x 6cm	4mm x 25mm	100 pcs/ Box
AFI6710	6.7 cm x 10cm (Hình oval)	7mm x 25mm	50 pcs/ Box